
REF PP17P

# PROPULSE® IRRIGATOR


UK  
CA  
0086

CE  
2797

**mirage**  
Health Group

2025-08-14 v 7.0

# PROPulse®

Propulse® ProScoop™  
INS0040

Propulse® Cleaning Tablets  
(x200 per pack)  
CL0001

Propulse® QrX™ Single  
Use Tip (x100 per pack)  
KIT5002

Propulse® ProTect™ Cape  
KIT1012

Propulse® Carry Case  
KIT1005

## INTENDED PURPOSE

The Propulse® Ear Irrigator is intended to:

- a) Facilitate the removal of cerumen and foreign bodies that are not hygroscopic from the outer aural meatus by irrigation with warm water.
- b) Remove discharge, keratin or debris from the external auditory meatus by irrigation with warm water.

Reasons for using this procedure are to:

- a) Correctly treat otitis externa where the meatus is obscured by debris.
- b) Improve conduction of sound into the ear, where impacted wax is believed to be the cause of a hearing defect.
- c) Examine the external auditory meatus and the tympanic membrane.
- d) Remove a cause of discomfort.

This procedure should ONLY be carried out by a suitably qualified healthcare professional.

## WARNINGS AND CAUTIONS

- This manual must be read and understood before the Propulse® Ear Irrigator is used. Failure to follow manual may lead to Device malfunction or ineffective ear irrigation.
- Only suitably trained staff should use the device. Mirage can advise on the availability of training courses offered by relevant organisations.
- The Propulse® QrX™ Tip is "Single Use" and should be disposed of in accordance with local authority guidelines after use.
- Reuse of the Propulse® QrX™ Tip increases the risk of cross-infection.
- Propulse® Ear Irrigator must not be immersed in water.
- Only clean the device as specified in this manual (See page 10).
- If any changes in performance occur, turn off the Propulse® Ear Irrigator, disconnect from the mains electricity supply and DO NOT use (See page 11).
- The device has no user serviceable parts (See page 13).
- Use recommended Propulse® accessories only.
- Do not use Propulse® accessories with other devices.
- If the device is to be used for domestic visits, it is strongly recommended that a Propulse® Carry Case is used to prevent damage and contamination.
- The Propulse® Ear Irrigator is not user repairable and should be returned to your Propulse® supplier or Mirage Health Group (UK customers only) for service and/or repair. **It is recommended that the Propulse® Ear Irrigator is serviced annually.**
- Only the footswitch provided can be connected to the Propulse® Ear Irrigator.
- The power supply must not be positioned such that it is difficult to disconnect from the supply mains.
- No modification of this equipment is allowed.
- It is unsafe to use any accessories (including power supplies), detachable parts or materials not specified in this device in the instructions for use.
- Caution must be taken while use in home environment due to pediatric risk of choking from small parts of device & strangulation of cables for baby or child.
- Propulse® unit must be cleaned prior to first use removing the yellow sticker inside the tank. Failure to remove sticker prior to use will cause obstruction in fluid circulation .


- Do not apply Propulse® QrX™ tip to ear if they are above or below ambient temperature.
- While using Propulse® unit - it can be unsafe to
  - use accessories, detachable parts and materials not described in the instructions for use.
  - interconnect this equipment with other equipment not described in the instructions for use.
  - modify the equipment.

### Illustrative Photograph of Irrigation

Note: clinical procedures relating to the use of ear irrigators are provided in Pg No.13

### Illustration of Device Connection

**Please note:** Damage caused to your Propulse® Ear Irrigator by the use of accessories, consumables or service agents not recommended by Mirage Health Group, will invalidate your warranty.


**DO NOT IRRIGATE** the ears if:

- Consent is not given and/or the patient is uncooperative.
- Previous complications occurred following this procedure.
- There is a history of a middle ear infection in the last six weeks.
- The patient has undergone ear surgery (apart from grommets that have extruded at least 18 months previously and the patient has been discharged from the ENT dept.)
- The patient has a perforation or there is a history of a mucous discharge in the last year.
- The patient has a cleft palate (repaired or not).
- In the presence of acute otitis externa; an oedematous ear canal combined with pain and tenderness of the pinna.
- If patient complains of pain/discomfort - STOP EAR IRRIGATION IMMEDIATELY.

*\*Please refer to Contraindications on page 12.*

**COMPONENT/PARTS**

- Reservoir Tank (detachable part)
- Lid of Reservoir Tank (detachable part)
- QrX™ Tip (accessory) (REF: KIT5666)
- Handle and Hose (non-Detachable part)
- Handle holder (non-detachable part)
- Waterflow/Pressure control switch
- On/Off switch
- Footswitch (accessory) (REF: KIT5016)
- Mains power adaptor (accessory) (REF: KIT5015)
- LED Indicator Light

The Propulse® Ear Irrigator consists of:

- The main unit and the following user controls:
  - An On/Off switch.
  - LED Indicator light - to indicate that the equipment is ready for normal use/ output.
  - A footswitch (REF: KIT5016) which (when pressed) starts the flow of water. The water stops when the footswitch is released.
  - A mains power adaptor. (REF : KIT5015)
- Water container/reservoir (1) is removable to facilitate filling and cleaning. Horizontal line indicates the correct level of water required for normal use, as well as the correct level of water required to dissolve cleaning tablet.
- Mushroom Valve (Detachable part) – to retain the water in the reservoir when it is removed from the Propulse® machine.
  - Handle and non-detachable hose. The Handle accommodates the Propulse® QrX™ Single Use Tips (REF: KIT5666).
- The footswitch - is connected to the main body via a jack plug/socket connection. The device will only function if the footswitch is connected.
- Instructions for use.

**Please note:** Residual water in the handle and hose will continue to flow if the handle is not held in the vertical position or, if the handle is held in a position that is lower than the machine. To prevent residual flow, it is recommended that the handle be returned to its holder on the machine.

## TECHNICAL DATA

Performance:

Description of function:

The device functions by pumping water from the reservoir through a hose, to the handpiece and the Propulse® QrX™ Single Use tip, using an electrically powered piston pump controlled by the footswitch. The device has a valve which adjusts the flow rate.

Flow rate:

Adjustable between 80 and 300ml/min

Applied part:

The Propulse® QrX™ Single Use tip

Applied part classification:

Type BF applied part

Physical characteristics:

### Irrigator

Mass: 625 g

Dimensions: 115mm X 220mm X 110mm

Reservoir capacity: 500ml

IP Rating: IP 21

### Footswitch (Supplied as an accessory)

Mass: 105g

Dimensions: 100mm X 60mm X 35mm

IP Rating: IP 21

### Power Adaptor

Mass: 135g

Dimensions: 75 x 85 x 48 mm

IP Rating: IP XO

Duty cycle:

Non-continuous use/operation:

Maximum of 10 minutes use followed by a minimum 2 hours of non-use

Power adaptor:

(Manufactured by ELijintek Inc)

This device is intended to be used with a separate power source, forming a Medical Electrical System.

Medical Power supply

Product Ref: KIT5015

Protection Class II

Input 100-240V ~ 50/60Hz

Max 0.45A

Output 9V DC2A

Note: In case of emergency, removal of the separate power supply mains plug from the supply mains is used as the means to isolate the ME system from the supply mains.

Continuous operating conditions

Temperature: 5 – 40°C

Relative Humidity: 10 – 90%

Pressure Range: 700 hPa to 1060hPa

Transport and storage between uses:

Temperature range : – 25 °C to + 5 °C, and

+ 5 °C to + 35 °C at a relative humidity up to 90 %, non-condensing;

> 35 °C to 70 °C at a water vapour pressure up to 50 hPa

Pressure range: 700 hPa to 1060 hPa

Humidity range: 10 % to 90 %

Equipment remains operational after transport or storage in the above environmental range and after having been removed from its protective packaging (Propulse Carry Case REF: KIT1005) and subsequently between uses.


Electrical Safety:

EN60601-1 and EN60601-1-11

EMC Compliance:

EN60601-1-2

## Symbols Used in IFU & Device.

	Caution		Direct current
	Type BF applied part		Variable flow rate
<b>IP21</b>	<b>IP- Ingress Protection Rating</b> 2-Protected against solid foreign objects of 12.5mmØ and greater. 1- Protected against vertically falling water drops .	10min/2h	Duty cycle
	Do not re-use		Catalogue number
	CE mark		Lot or batch number
	UKCA mark		Indoor use only
	To indicate that separate collection for waste electric and electronic equipment (WEEE) is required		Indicates the item is a medical device
	Manufacturer		Authorised representative in the European Union
	Power on		The instructions for use must be read
	Power off		Operator's manual; operating instructions
	Permitted humidity limits		Permitted temperature limits
	Date of manufacture		Unique Device Identifier
	Foot Control (Foot switch )		Serial number
	Polarity of d.c. power connector*		Atmospheric pressure limitation
	To indicate the number of pieces in the package		Class II equipment
<b>PROPULSE</b>		Mirage Health Group, Propulse Registered traded logo mark for traceability and part/accessory identification.	

## OPERATING INSTRUCTIONS

These instructions are for general use. When required, refer to the detailed information in the second half of this manual.

- Ensure that ONLY suitably trained clinicians operate the device.
- Ensure that the warnings and cautions are observed.
- Ensure that the patient exhibits no contraindications (refer to page 12).
- Ensure the unit has been cleaned prior to first use, and every day prior to use (refer to page 10 for detailed cleaning guidance).
- The device may only be operated whilst connected to the mains electrical supply using the Power Adaptor and footswitch provided. These connect to the sockets on the side of the main unit. Ensure that the footswitch and power lead do not cause a trip hazard
- The reservoir should be removed prior to filling.
- The water reservoir must be filled to the horizontal line on the front of the reservoir with warm water. Note water which is too warm or too cold may cause vertigo as a side effect.
- Monitor temperature regularly to ensure that patient comfort and safety is maintained. Refill as necessary.
- Fit a new Propulse® QrX™ Single Use Tip to the Handle.
- The flow of water is controlled by depressing the footswitch. The flow stops when the switch is released.
- Position the patient in an upright chair, close to the worksurface the Propulse Ear Irrigator is placed on allowing ease of access to the ear. It is recommended that a Noots tank (Not Supplied) is used to collect water draining from the ear and while testing the flow and temperature of water.
- Turn the On/Off Switch (7) to the position marked "I".
- Adjust the Waterflow Switch (6) to an appropriate level.
- Adjust the waterflow to an appropriate value that is compatible with the treatment requirements and patient comfort.
- Direct the irrigator tip into the Noots tank or similar receptacle and depress the footswitch for 10-20 seconds in order to circulate the water through the system and eliminate any trapped air or cold water.
  - the time from switching "on" until the me equipment is ready for normal use: 10-20 seconds.
  - the time required for the equipment to warm from the minimum storage temperature between uses until the equipment is ready for its intended use when the ambient temperature is 20 °C: 10-20 seconds.
  - the time required for the equipment to cool from the maximum storage temperature between uses until the equipment is ready for its intended use when the ambient temperature is 20 °C: 10-20 seconds.
- Repeat the procedure until the water jet reaches the desired clinical temperature.
- This procedure should also be repeated for 5 seconds every time a new QrX™ Single Use Tip is used.
- Ensure the water is warm before presentation to the patient.
- During treatment you can pause the flow by releasing the footswitch.
- After treatment empty the reservoir and operate the device to purge any residual water.

- Residual water, and water collected in the noots tank can be disposed of with normal waste water.
- Remove the Propulse® QrX™ Tip and dispose of in accordance with local authority guidelines.
- Turn off the On/Off switch after use **and** disconnect from the power supply.
- Clean the Propulse® Ear Irrigator unit **every** morning **prior** to use, with a Propulse cleaning tablet (See page 10).
- The Propulse® Ear Irrigator should only be transported in a Propulse® approved carry case to prevent damage or contamination.
- If any changes in performance occur, turn off the Propulse® Ear Irrigator, disconnect from the mains electricity supply and **DO NOT** use. Please refer the machine to Mirage.
- Do not touch the power supply output connector and patient simultaneously
- Do not touch the footswitch connection port and patient simultaneously

## DETAILED INFORMATION

### Fitting the Footswitch

The Propulse® footswitch is connected to the main unit by a socket in the side of the device. The Propulse® Ear Irrigator will NOT operate unless the footswitch is connected.

### Fitting the water reservoir and filling the water

It is recommended that:

- i) The water reservoir is removed from the device for filling and that the lid is always in place when the water container is in-situ on the device.
- ii) The water reservoir should be filled to the horizontal line on the front. This helps to eliminate the risk of spillage.
- iii) Water at 38-40°C is recommended. Higher temperatures increase the risk of burns to the patient. Lower temperatures increase the risk of patient discomfort. Both high and low temperatures can cause vertigo.
- iv) While fitting back water reservoir tank and lid ensure there is no visible leakage of water before equipment is switched on /operation.

### Fitting the Propulse® QrX™ Tip

The Propulse® Ear Irrigator is designed to be used only with Propulse® QrX™ Single Use Tips. Use one Propulse® QrX™ Tip per treatment.

### To fit a Propulse® QrX™ Tip

- 1) Remove Tip from packaging - Tips are non-sterile.
- 2) Push Tip into Handle until a click can be heard.

### To remove a Propulse® QrX™ Tip

- 1) Retract the QrX™ Locking Collar using a thumb.
- 2) Grasp the used QrX™ Tip between forefinger and thumb and pull gently from the QrX™ Handle.
- 3) Dispose of the used Tip in accordance with local authority guidelines. DO NOT REUSE TIPS.

Propulse® QrX™ Tips are available to purchase in boxes of 100 individually wrapped (non-sterile) Tips from your normal Propulse® supplier or from Mirage directly (UK only). Propulse® QrX™ Tips are clearly branded with the Propulse® logo on the Tip and its packaging. Only branded Propulse® QrX™ Tips should be used with the Propulse® Ear Irrigator.

### Replacing the Mushroom Valve

The Mushroom Valve is specifically designed to prevent water flowing out of the reservoir whilst filling. Should a replacement Mushroom Valve be required, please follow the steps below -

**Firstly, ensure you are fitting the correct mushroom valve for this model Propulse® Ear Irrigator.**

This will help prevent possible damage to the water inlet valve.

- 1) Remove the reservoir from the Propulse® Ear Irrigator.
- 2) Remove the old Mushroom Valve from the reservoir and discard
- 3) Insert a new Mushroom Valve without bending or exerting undue force on the legs of the Mushroom Valve, into the reservoir.
- 4) Check the condition of the 'O' ring on the base of the reservoir and if worn replace with a new unit specific 'O' ring.
- 5) Return the reservoir to the machine

## Cleaning Instructions

The importance of using the correct strength cleaning solution cannot be overstated. A solution that is too strong will in time damage the Propulse® Ear Irrigator. A solution that is too weak will fail to provide the correct level of cleaning. Mirage Health Group recommends the use of Propulse® CHLOR-CLEAN® Tablets. They are easy and effective to use and provide a measured / specific strength of cleaning solution that is safe and kind to the Propulse® internal components.

### **Ensure the unit has been cleaned prior to first use.**

1. Place warm tap water into the reservoir up to the horizontal line on the front.
2. Place one Propulse® CHLOR-CLEAN® tablet into the reservoir and allow it to dissolve completely.
3. Once dissolved run the machine until the cleaning solution leaves the handle. This ensures that the cleaning solution has reached all of the internal components.
4. Leave the solution in place for 10 minutes.
5. After 10 minutes remove the reservoir with the remaining cleaning solution and discard.
6. Fill the reservoir with clean, well run, cold tap water and return to the Propulse®.
7. Run the Propulse® ensuring that all remaining cleaning solution has been flushed through.
8. Remove reservoir, discard water and dry reservoir thoroughly with a lint free cloth.
9. Return reservoir to the Propulse® - it is now ready for use.

Do not attempt to clean the Propulse® QrX™ Tip. Use one Propulse® QrX™ Tip per patient treatment and discard to clinical waste after use as this reduces the risk of cross infection between patients. Before using a new Propulse® QrX™ Tip on a patient the machine should be turned on for at least 5 seconds, as set out in the operating instructions.

External cleaning of the Propulse® Ear Irrigator should be undertaken by hand, wiping with a damp lint free cloth only. Apply liquids to the lint free cloth not the unit. Do not immerse the unit in water. Mild detergents and disinfectants may be used externally, do NOT use flammable, solvent or abrasive cleaning materials as these will damage the device. There is no need to clean the device between uses on multiple patients.

## **Power Adaptor**

Connect the outlet lead of the Power Adaptor to the Power Adaptor socket marked on the end of the product and to the mains electrical supply. Make sure the cord and Power Adaptor are positioned so that they will not be subjected to damage or stress or present a trip hazard.

Only use a Propulse® branded Power Adaptor.

To reduce the risk of electric shock, unplug the unit from the power source before attempting to clean it externally.

The Power Adaptor must not be used outdoors or in damp areas.

Removal of the separate power supply mains plug from the supply mains is used as the means to isolate the ME system from the supply mains

The Propulse® Power Adaptor will have been supplied with an appropriate plug for your region OR a selection of international plugs. Please fit the appropriate plug for your region.

If there are any problems connecting to the mains electrical supply consult a qualified electrician.

## **Electromagnetic Disturbances**

The Propulse and its accessories have been designed and tested for compliance to EN 60601-1-2 for immunity and emissions of electromagnetic interference.

In the unlikely event of disturbance of other devices caused by the Propulse ensure that the wires of the Propulse are adequately separated from those of other devices and move the Propulse away from the other devices. If the problem persists contact the service department

As the Propulse contains minimal electronic circuitry, interference with the Propulse by other devices is extremely unlikely but if interference does occur follow the guidance as per the previous paragraph with respect to positioning of wires and devices.

## Contraindication to irrigation

## Rationale

If the patient has experienced any complications from a previous episode of irrigation with water.

If the patient did not tolerate a previous episode of irrigation it would be unwise to repeat the procedure in case the symptoms are exacerbated.

There has been evidence of a middle ear infection (otitis Media) in the last 2 months.

The tympanic membrane may be vulnerable to damage due to the adverse effect infected fluid may have on the ear drum.

If the patient has a broken skin with blood evident in the external ear canal

Patients existing broken skin in the EAC increases the risk of further bleeding and/or infection with irrigation

The patient has undergone any form of ear surgery apart from grommets, which are documented to be extruded from the tympanic membrane for over 2 years and the patient is discharged from the ENT department.

There will be a weakness to the structure of the ear canal and tympanic membrane after surgery. This does not include cosmetic surgery to the pinna (for example repair of bat ears). If the tympanic membrane is intact 2 years post grommet extrusion, there should not be an increased risk of damage to the tympanic membrane.

There is a suspected or actual perforation present or there is a history of mucous discharge from the ear in the last 2 years.

A mucous discharge would indicate a perforation and water entry under pressure could cause infection or damage the delicate middle ear structures.

If the patient has a cleft palate (regardless of whether it has been repaired or not).

A cleft palate indicates an underdeveloped facial skeleton and as such the tympanic membrane and middle ear structures could be more vulnerable to damage.

In the presence of acute otitis externa (pain, swollen ear canal and tenderness of the pinna).

Although it is essential to thoroughly clean the infected ear canal, when it is swollen debris should be removed by microsuction.

Profound hearing loss in one ear.

There is a risk associated with any intervention and when a patient completely relies on one ear for hearing (as the other ear has a profound hearing loss) any risk to this ear is unacceptable.

## Caution when irrigating with water in the following groups of patients

## Rationale

Patient is taking anti-coagulants.

The lining of the ear canal is delicate and there is a higher risk of bleeding so ensure trauma to the ear canal is avoided.

The patient is diabetic.

The pH of wax in patients who are diabetic is a higher pH than average, increasing their vulnerability to infection.

Tinnitus.

Although wax impaction can cause tinnitus, trauma to the tympanic membrane may exacerbate this.

Vertigo.

This is also a symptom of wax impaction but irrigation can trigger an episode so ensure appropriate water

Radiotherapy that has involved the ear canal.

A radiated ear canal can develop bony necrosis so wax should be removed before it becomes hard and trauma to the canal should be avoided.

## Maintenance & Safety Inspections

**To ensure optimum performance the Propulse® Ear Irrigator should be serviced every 12 months. Service or repairs conducted by unauthorized agencies/organisations are not permitted and invalidate any or implied warranties from Mirage.**

The Propulse® Ear Irrigator should under-go routine electrical safety testing to ensure that it remains safe to use, in accordance with EN ISO 62353:2014.

Users of the Propulse® Ear Irrigator should carry out regular inspections to ensure that the handle and hose, power adaptor and cable, reservoir, footswitch and main body of the machine are free from damage prior to use. If any damage is evident, the Propulse® Ear Irrigator should NOT be used until replacement parts have been fitted. These parts have the same service life as the rest of the device.

Only Propulse® branded items should be used with the Propulse® Ear Irrigator.

The Propulse® Ear Irrigator is not user repairable and should be returned to your Propulse® supplier or Mirage Health Group (UK customers only) for service and/or repair:

### Mirage Health Group

11 Tewin Court, Welwyn Garden City, Hertfordshire  
AL7 1AU  
UK  
Tel - +44 (0) 845 130 5445

No part of the device can be serviced or maintained while in use with a patient.

All parts of the device are suitable for use within the patient environment.

The expected service life of this device and its accessories (excluding the single-use Propulse® QrX™ Tip) is 4 years.

The expected shelf life of single-use Propulse® QrX™ Tip is 10 Years

The clinical procedures relating to the use of ear irrigators can be found on the following website:

[www.earcarecentre.com](http://www.earcarecentre.com)

[www.entnursing.com/earcare.htm](http://www.entnursing.com/earcare.htm)

Mirage is not responsible for the content or maintenance of third-party internet sites.

Mirage can also advise on the availability of training courses offered by the relevant organisations.

Additional information on use of the Propulse® can be found at:

<http://www.youtube.com/user/MirageHealthGroup>

## Warranty

The Propulse® Ear Irrigator carries a twelve-month warranty (\*subject to conditions) from the date of original purchase. Should any defect arise due to faulty material or workmanship, Mirage Health Group will, upon receipt of the faulty Propulse® Ear Irrigator, proof of purchase, information relating to the nature of the fault and details of where the item was purchased, rectify the fault at no cost to you.

Should any of the “Accessory” items (listed below) prove to be faulty as a result of defective material or workmanship, Mirage Health Group will rectify the issue free of charge upon receipt of the faulty accessory (\*subject to conditions).

“Accessory” items are: Footswitch; Reservoir / Tank and Lid; Mushroom Valve and Washer; QrX™ Tip.

\*Conditions of Warranty (applicable to Propulse® Electronic Ear Irrigator and “Accessory” items).

### The warranty does not cover:

- Accidental damage or damage caused by misuse.
- Faults caused due to lack of maintenance.
- Damage caused by using the Propulse® Ear Irrigator for any use other than its intended use.
- Damage caused as a result of repair by any unauthorized agents - ONLY Mirage Health Group should undertake repairs.
- Damage caused by the use of accessories / cleaning products that have not been recommended by Mirage Health Group as suitable for your model irrigator.

This warranty is in addition to, and does not diminish your statutory or legal rights.

Additional user manuals and other accessories are available from Mirage Health Group Ltd at:

Mirage Health Group  
11 Tewin Court, Welwyn Garden City,  
Hertfordshire AL7 1AU UK

Tel - +44 (0) 845 130 5440

Fax - +44 (0) 845 130 6440

[www.miragehealthgroup.com](http://www.miragehealthgroup.com)

[uksales@miragehealthgroup.com](mailto:uksales@miragehealthgroup.com)

[internationalsales@miragehealthgroup.com](mailto:internationalsales@miragehealthgroup.com)

EC	REP
----	-----

**Medical Device Management Ltd**  
Block B, The Crescent Building,  
Northwood, Santry,  
Dublin 9, D09 C6X8,  
Ireland

## **Environmental Protection**

This symbol on the products and/or accompanying documents means that used electrical and electronic products should not be mixed with general waste. Please return to Mirage Health Group or dispose of through locally approved disposal service for electronic equipment. Penalties may be applicable for incorrect disposal of this waste, in accordance with national legislation.

### **Disposing of used Propulse® QrX™ Tips**

Disposal should be in accordance with local authority guidelines and regulations for the disposal of clinical waste. Propulse® QrX™ Tips should not be disposed of in municipal waste.

## **Transportation**

Before the Propulse® Ear Irrigator is transported, the reservoir must be emptied and the machine should be operated until the handle and hose are empty of liquid. The reservoir should then be dried using a paper towel.

For safe transportation of Propulse® Ear Irrigator, Mirage recommends that the Propulse® Carry Case is used to prevent damage or contamination.

For internal cleaning, please refer to page 10.

## **Serious Incidents**


If any serious incident has happened where this device is suspected to be a contributory factor, then this should be reported to Mirage Health Group, using the contact details listed on the previous page, and the competent authority in the country the device was used in.

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Extra gebruikershandleidingen en accessoires zijn beschikbaar van Mirage gezondheid Group Ltd op: